# Femcon® Plus

Vaginalkonen

Vaginal Cones

Cônes vaginaux

Coni vaginali

Conos vaginales

Cone vaginal

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**Please study the instructions carefully, for you will only achieve the results desired if you use Femcon® Plus properly.**

### **Incontinence – a problem you can solve**

One woman in three suffers from incontinence. But this „muscle weakness“ is not only a physical problem. As the uncontrollable loss of urine occurs in particular when the person affected laughs, sneezes, coughs, runs or puts some strain on her body, many women can no longer enjoy their lives the way they used to.

### **Reasons for incontinence**

In most cases, incontinence was triggered off by having given birth to one or more children. Research shows that nerves connected with the pelvic muscles are often harmed by giving birth. In the course of time, this can lead to an increased weakness of all pelvic muscles. The muscles then are no longer able to close the bladder against the urethra. The female pelvic muscles may also lose their flexibility because of weak connective tissue, hard physical labour, overweight or hormonal changes during menopause.

*It is important that you exercise daily and over a long period of time.*

### **Preventive steps**

In order to prevent incontinence caused by stress, it is very important to exercise daily to develop the pelvic muscles – and you should start in the weeks following childbirth, of course under a specialist's supervision. After the lochia has subsided, well-directed exercises with Femcon® Plus Vaginal Cones are the best way to strengthen the strained pelvic muscles. In the beginning, it is enough to exercise 8 to 10 minutes twice a day. Ageing women or women who have just given birth can often feel that their vaginas don't contract properly any more. Well-directed exercises with Femcon® Plus help rebuild active, firm muscles. To make the use of Femcon® Plus Vaginal Cones as comfortable as possible and to enhance their effectiveness, we recommend that you apply Ladysoft gel.

## **Steps to follow if you suffer from incontinence caused by stress**

Only pelvic muscles that work properly are able to keep the organs of the small pelvic (bladder, uterus) in the right position. By exercising carefully and continuously you can get rid of, or at least improve considerably, any form of incontinence. Surgery can be avoided in almost 70 % of all cases of incontinence (slight to moderately severe) and first/second degree prolapse of organs.

## **How Femcon® Plus Vaginal Cones work**

After having put some Ladysoft gel on the cone, insert it into your vagina as if it were a tampon. Its weight and smooth surface convey a „slippery“ feeling and you think it will slide out. This feeling triggers off a nervous reaction. The pelvic muscles contract around the cone, stop it from sliding out and push it back up („biofeedback phenomenon“). By tightening the muscles intentionally, they can be rebuilt. Thereby you can test your pelvic muscles and strengthen them by well-directed exercises.

## **Simple usage, great results**

How to use Femcon® Plus Vaginal Cones Before you start exercising with Femcon Plus Vaginal Cones you should talk to your doctor or physiotherapist about how the pelvic muscles work. If you still cannot tighten your muscles after having tried several times, we advise you to talk to your gynaecologist. He will know how to make you aware of these muscles.

Exercising with Femcon® Plus Vaginal Cones is easy, and if you use them continuously you will feel the results after a few weeks. Femcon® Plus Vaginal Cones come in sets of five. The cones are shaped like tampons and have the same size but vary in weight. A tube of Ladysoft gel completes the package. The cones have numbers that tell you the weight (1/20 g, 2/32,5 g...5/70 g).


Don't forget to wash your hands before you start the exercise, then apply two small drops of Ladysoft gel on the cone or at the entrance of your vagina. Thus you can insert the cone more easily and enhance the effectiveness of the exercise. Insert the cone with the least weight (cone no.1) deep into your vagina – you should be able to hold it there at least a minute without any effort (while standing upright). If you have completed this exercise successfully, take cone no.2 and repeat the procedure, then cone no.3 and so on. The cone that you can just about hold informs you about the condition and strength of your pelvic muscles.

*Slowly but steadily on the road to success.*


Always start exercising with the heaviest cone you were able to hold without any effort the last time. Exercise several times every day while standing or walking about. During this time, you don't have to suppress the need to cough, sneeze or laugh or any other minor workouts. Initially, you might not meet your own demands, but don't be discouraged. Insert the cone again, and concentrate hard on tightening your pelvic muscles. If you're able to hold the cone for at least 8 to 10 minutes several days in a row, you repeat the exercise with a heavier cone in the same way.

## Femcon® Plus – easy to use

1. Put some drops of Ladysoft gel on the cone and insert it gently


2. If your pelvic muscles are weak the cone slides downward

3. Exercises and their effects


A

Because of its weight, the cones slides downward


B

If you tighten your pelvic muscles...

C

...the cone is pushed back up

## **A successful treatment**

Several clinical studies have proven that daily and careful exercises reduce incontinence distinctly and strengthen the pelvic muscles effectively.

After a few weeks, with the help of Femcon® Plus Vaginal Cones, you learn to train your pelvic muscles even without using a cone. If there's a significant improvement after having practised for 6 to 8 weeks, you can interrupt your exercise programme. You can check whether your muscles work properly from time to time with Femcon® Plus Vaginal Cones and, if necessary, you simply exercise intensely about 2 to 4 weeks. With this kind of therapy the problem won't reoccur in case your connective tissue shows any signs of weakness. A treatment with Femcon® Plus Vaginal Cones can be carried out successfully at any age.

## **Comments**

During menstruation or any inflammation of the bladder or vagina it is advisable to interrupt the exercises with Femcon® Plus Vaginal Cones. When you start exercising, a feeling of stiffness in the area of the pelvis and intestines can occur. This is just a sign that the muscles start building up.

It is very important to clean used cones regularly with soap and warm water. To disinfect the cones, just use a regular disinfectant available on the market. Before you reuse a cone, be sure to rinse off all the disinfectant. The cones can be sterilized with gas or liquids, but not with steam or hot air (that would deform the cones).

Femcon® Plus Vaginal Cones are not to be used on children. Teenagers are strongly advised to consult a gynaecologist. It is not recommendable to use Femcon® Plus Vaginal Cones and we advise you to consult your gynaecologist if your vagina is very short (e.g. in case of atrophy caused by age, craurosis fornicis vaginae or surgery) or in case of a prolapse of the vagina or other organs.

## Femcon® Plus Vaginal Cones

Femcon® Plus Vaginal Cones are available at Pharmacies as well as in special medical stores or at:

medesign I.C. GmbH  
Dietramszeller Str. 6 · D-83623 Dietramszell  
[www.medesign.de](http://www.medesign.de) · [info@medesign.de](mailto:info@medesign.de)  
Phone +49 (0)8027 90738-0 · Fax +49 (0)8027 90738-27

Usable until

Item number (MED...)

LOT

Attention

Date of construction

Read and follow the instructions for use

Producer

Single patient – multiple use

0483 CE-mark with identification number of the notified body

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medesign I.C. GmbH · Dietramszeller Str. 6 · D-83623 Dietramszell